

# PHARMACEUTICAL DOSAGE FORMS

Tablets

In Three Volumes

VOLUME 1

EDITED BY

Herbert A. Lieberman  
*Warren-Landell Company, Inc.  
Johans Plains, New Jersey*

Leon Lachman  
*United Laboratories, Inc.  
Makati, Metro Manila  
The Philippines*

MARCEL DILKER, INC. New York and Basel

MARSHALL 10971



MARSHALL 11189

212

Figure 11




Figure 11. Schematic of a layer press. Refer to text. (Thomas Engineering.)

# Compression-Coated and Layer Tablets

When a layer is ejected, the upper tamping roll is lowered slightly to exert more pressure upon the layer. This action will prevent damage to the layer as it strikes the take-off blade and is directed into the collection box. Once the layer strikes the take-off blade and is directed into the collection box, the upper punch strikes the next filling station, they are quickly pulled down by a lowering ram so that they are not struck by the upper punches. The latter are always descending into the dies to make the next tamping or compression stroke.

The leading and trailing edges of each feed frame are equipped with wipe-off blades which direct any powder that escapes from the feeders into collection boxes. The blade on the trailing edge of the first feed frame guides the completed tablets away from the die into the collection bin. Vacuum tubes at each filling unit suck away any powder or granulation that remains in the lower punch lines during weight checks. Although the punches are raised flush with the die table at this time and do not drop as they pass under the feed frames, they do trap a small amount of granulation in the depressions in their tips.

If an adjustment in the weight or thickness of the first or second layer is necessary, then the weight of each succeeding layer will probably need correction, since weight is related to the fill volume.

The second type of machine is similar to the one described above, except for the manner in which weight checking is handled. Instead of a ram arrangement for ejecting the layers, the pressure on the first layer is increased, and the layer is made so hard that the next layer will not bond to it. Thus both layers are easily separated for weighing. This effect is achieved by actuating a pneumatic cylinder which raises the lower tamping roll. There is an adjustment to control the distance which raises the lower tamping roll. There is an adjustment to control the distance that the compression roll may rise. Engraved or engraved upper punches provide a key for each layer and need to hold them together. Gentle shaking may be required to separate the layers in this case. Table 3 provides specifications for several typical layer presses currently available.

## V. Formulations (Layers)

As with compression-coated tablets, the granulation for layer tablets should be readily compressible for good bonding between layers. Doublet dies should be kept in a minimum; the less dust, the cleaner the scrape-off at each feed frame. It may be necessary to separate out that fraction of a granulation which is finer than 70 or 80 mesh. Such material is not discarded but added to the next lot and regranulated. Just enough, however, must be finely divided, their efficiency determined by the degree of fineness. Since the lubricant must cannot be avoided, the quantities used should be kept minimal. The metallic stearate present on the additional difficulty is that they interfere with the bonding of the layers. Stearic acid and the hydrogenated fats are better lubricants from this point of view. Nevertheless, granules should be small, less than half the thickness of the layers; otherwise, the lines of demarcation between layers will be uneven.

Equal weights of granulation will not necessarily lead to equal thickness of the layers. That will depend on the compression ratio of the formulations. It is easy to compensate for by adjusting the weights required for each layer. (It is not necessary, however, that each layer have the same thickness.) The shape of the punches also plays a role: punches with beveled edges or concave faces will make the top and bottom layers of a three-layer tablet appear thinner than the middle one. Flat-faced tooling will produce equal thickness of the layers, but

MARSHALL 11190

MARSHALL 11191

218

Table 2  
Specifications for Snow Layer Process

Specification	Manufacturers and model designation				
	Manufacturer 39	Thickness 47	Pressure 65	Quantity 61	Notes 805-2
Number of dies	39	47	55	61	29 45
Maximum pressure (lb/in. <sup>2</sup> )	0.1/2	0.1/2	0.1/2	4	8.1/2 4
Maximum tablet diameter	5/8 in.	5/8 in.	5/8 in.	7/16 in.	3/8 in. 5/8 in.
Depth of fill	11/16 in.	11/16 in.	11/16 in.	11/16 in.	15 mm 11/16 in.
Maximum layer thickness (prior to pressing)	1/4 in.	1/4 in.	7/16 in.	7/16 in.	7 mm 1/2.0 in.
First layer	1/4 in.	1/4 in.	1/4 in.	1/4 in.	7 mm 1/4 in.
Second layer	1/4 in.	1/4 in.	1/4 in.	1/4 in.	7 mm 1/4 in.
Third layer	1/4 in.	1/4 in.	1/4 in.	1/4 in.	7 mm 1/4 in.
Maximum weight (tablets per minute)	1250	1500	5000	5000	417 2400

# Compression Coated and Layer Tablets

219 70a



Figure 12. Cross sections of layer tablets.

Unfortunately the edges of the tablets will tend to chip readily. Figure 13 shows cross sections of layer tablets and illustrates how the shape of the upper punch determines the shape of the layers. If the upper punch faces have monograms or other markings, the bonding between layers will be strengthened because the devices will act as keys between the layers. Additionally, precompression lengthens dwell time and aids in bonding. The formulae previously given for compression-coated tablets will serve as a guide for the development of formulations for layer tablets, with the exception of two of those for direct compression (Examples 1 and 2), which are composed entirely of fine substances.

An illustrative formula is one for an analgesic-antipyretic demulgent containing aspirin and phenylpropanolamine. A thin layer of placebo is placed between them to negate the chemical incompatibility of the active ingredients.

## Example 1: First Layer of Analgesic-Antipyretic Demulgent

Ingredient	Quantity per Tablet
Phenylpropanolamine HCl NF	15.16 mg
Lactose USP	76.50 mg
Sucrose HSP	2.00 mg
Vale USP	2.00 mg
Tragacanth USP	2.00 mg
Polyethylene glycol 6000 USP	4.00 mg
Purified water USP	q. s.
Anhydrous alcohol	q. s.

Screen where necessary to break down agglomerates or lumps (20 mesh screen is satisfactory) and blend the phenylpropanolamine, lactose, sucrose, vale, and tragacanth. Combine the purified water and alcohol and dissolve the polyethylene glycol in the mixture, using heat (50 to 60°C). Add this solution to the mixed powders. Continue mixing until the mass is evenly moistened and granular. Use additional water-alcohol, if necessary.

Dry in an oven with circulating, dehumidified air at 40°C. Pass the dried granules through a 20 mesh screen on a Fiamillo or Tornado mill running at medium speed with knives forward.

MARSHALL 11192

Example 16: Second Layer of Analgesic-Antipyretic  
Decongestant

Ingredient	Quantity per tablet
Lactose USP	48,000 mg
Confectioners sugar USP	24,000 mg
Starch USP	7,070 mg
TSPC celite	9,095 mg
Purified water USP	q.s.
Sodium acid USP	0,995 mg

Pass the lactose and confectioners sugar through a 20 mesh screen and blend them with the starch in a suitable mixer. Dissolve the color in the water and add to the mixed powders. Continue mixing until the color is uniformly dispersed and the mass is granular. Dry in an oven at 40 to 45°C to a moisture content of 5% or less. Pass the dried material through a 20 mesh screen on a Torndahl comminuting machine. Return the granules to the mixer and add the stearic acid. Mix for 10 min.

Example 17: Third Layer of Analgesic-Antipyretic  
Decongestant

Ingredient	Quantity per tablet
Aspirin-Starch (20 mesh granules, 10% starch)	80.0 mg
Talc USP	10.0 mg

Blend in a suitable mixer until homogeneous (10 to 15 min).

Compress the three layers together using 3/8-in. diameter, flat-faced, beveled-edge punch. The weight of each layer is:

First layer, 160 mg  
Second layer, 80 mg  
Third layer, 100 mg

The top layer is the last layer to be pressed. Since it is the aspirin portion, it will be most resistant to extrusion from the dies.

Layer presses find employment in the manufacture of chewable antacid tablets. A possible formula for such a product follows. The material provides pleasant mouth-feel and sweetness, and the saccharine masks the latter. Peppermint flavoring has a long and honorable association with antacid preparations. The

MARSHALL 11193

# Cooperation Coated and Layer Tablets

221

sucrose units as the binder, although, of course, it also contributes to the taste of the tablets.

## Example 15: First Layer of Chewable Almond Tablet

Ingredient	Quantity
Magnesium oxide, heavy	200.0 mg
Mannitol NF	400.0 mg
Sucrose USP	60.0 mg
Saccharin sodium	1.0 mg
Purified water USP	q.s.
Magnesium stearate USP	7.0 mg
Peppermint oil	4.0 mg

Blend the magnesium oxide, mannitol, and saccharin in a double-arm mixer. Dissolve the sucrose in double its weight of water and add to the blended powders. Continue mixing until a moist, granular mass is formed, using additional purified water if necessary. Pass the batch through a 20 per-  
forated plate on a Fitzmill operating at low speed with beaters forward. Spread the material on trays and dry at 60°C. Pass the dried granules through a 12-mesh screen on a Fitzmill running at medium speed with knives forward. Return the granules to the mixing machine and add the peppermint oil. When the oil has been thoroughly dispersed, add the magnesium stearate. (If the oil is not added before the lubricant, the tablet will have oil spots on its surface.) Compress the layer at 677 mg using 5/8-in. diameter punches with flat tops and beveled edges.

## Example 16: Second Layer of Chewable Almond Tablet

Ingredient	Quantity
Aluminum hydroxide (dried gel)	200.0 mg
Mannitol USP	400.0 mg
Saccharin sodium	0.6 mg
Starch USP	22.0 mg
Purified water	100.0 mg
Oil of peppermint	3.4 mg
Magnesium stearate USP	7.0 mg
Color	q.s.

MARSHALL 11194



Blend the aluminum hydroxide, mannitol, and saccharin. Dissolve the color in the water and add the starch. Heat the mixture on a waterbath until the starch jells and forms a paste. Use the paste to granulate the blended powders. Add more water, if necessary, to form a lumpy mass. Pass the mass through a #8 perforated plate on a Fildmill running at low speed with the hammers forward. Spread the material on trays and dry at 55°C. Pass the dried granules through a 12 mesh screen on a Fildmill running at medium speed with knives forward. Return the granules to the mixer. Add the flavor first and then the magnesium stearate. Compress at 540 mg onto the first layer.

From the patent literature [14] there is this example of a three-layer tablet.

Example 20: Bottom Layer of Three-Layer Tablet

Ingredient	Quantity
Aspirin/Aspirin acid	210.0 g
PDC Yellow No. 5	4.0 g
Cornstarch	39.0 g
Talc	10.0 g
Chloroform	q.s.

Mix thoroughly and pass the mixture through a hammer mill. Add sufficient chloroform to obtain a wet granulation. Reduce the granules to a range of 20 to 40 mesh and dry overnight at a temperature of 120 to 140°F.

Example 21: Middle Layer of Three-Layer Tablet

Ingredient	Quantity
Phenacetin	100.0 g
Caffeine	15.0 g
Phenyltoloxamine dihydrogen citrate	10.0 g
Cornstarch	4.0 g
Powdered sugar	0.4 g
Distilled water	5.3 g
Magnesium stearate	3.0 g

Blend the phenacetin, caffeine, and phenyltoloxamine dihydrogen citrate. Prepare a paste by heating the starch and sugar in the water. Add the paste to the powders and form granules. Dry the moist mass overnight at 120 to

MARSHALL 11195

140° F. Reduce the mass to granules of about 20 mesh. Blend the granules with the magnesium stearate.

#### Example 32: Top Layer of Three-Layer Tablet

Ingredient	Quantity
Potassium phenethylmalonate	175.00 g
FD&C Red No. 3	0.03 g
Chloroform	q.s.
Magnesium stearate	3.00 g

Blend the first two ingredients and pass them through a hammer mill. Add sufficient chloroform to make a hard rubber-like mass. Break up the mass and dry overnight at 130 to 140° F. Reduce the dried material to about 20 mesh granules. Blend the granules with the magnesium stearate.

Using a three-layer press, compress the bottom layer at 254 mg, the middle layer at 197.6 mg, and the top layer at 175.00 mg.

Today, FD&C Yellow No. 5 would not be used with acetylsalicylic acid because of the possibility of allergic reactions.

Although compression-coated and layer tablets are a modest fraction of solid oral dosage forms, they provide two additional alternatives in solving formulation problems. They tend to be more expensive to manufacture than other tablets (except tablet retextures) because of the multiple granulations needed and the slowness of the special process used.

#### References

1. Noyce, P. J., British Patent 859955 (1959).
2. Stokes, E. J., U. S. Patent 3,048,591 (1963).
3. DeLong Gum Company, British Patent 497,556 (1930).
4. Kilian, F., British Patent 194,901 (1907).
5. Cooper, J., Farnside, D., and Windhauser, J., U. S. Patent 2,827,313 (1958).
6. Wolff, J., U. S. Patent 2,707,124 (1955).
7. Blumhage, F., Zepanas, J., and Quack, M., J. Amer. Pharm. Assoc. (Sci. Ed.), 47,13:697-670 (1958).
8. Windhauser, J., and Cooper, J., J. Amer. Pharm. Assoc. (Sci. Ed.), 45,4:543 (1955).
9. Leuchman, L., Spicer, P., and Sylvestermatics, H., J. Pharm. Sci., 52,4:379-382 (1963).
10. Bonnell, C., U. S. Patent 3,048,586 (1962).
11. Buchwalter, F., Crenshaw, A., and DeMoro, M., U. S. Patent 3,321,094 (1964).

MARSHALL 11196

## Suggested Reading

Remington's Practice of Pharmacy, 1901 ed., Mack Pub., Easton, Pa., 1973.  
Rifenshel, W. A., Die Tablette. Edeka Cantor KG, Aulendorf i. Wuerth., Germany,  
1966.

MARSHALL 11197